

## SPORLASTIC TRIAC II

### *Dynamic Skoliosis-Orthotics in an open frame construction*

#### INDICATIONS

Idiopathic Scoliosis up to 35° Cobb

#### CONTRA-INDICATIONS:

- Scoliosis' contour higher than TH 7
- Scoliosis' contour TH 11/12, L 4/5

Paralysed-Scoliosis

#### MODE OF ACTION

Permanent growth manipulative correction pressure on the spine while lying, sitting and standing  
 Stabilization of lateral and rotational deviations

Minimal imposition – in physical (skin) and psychological (aesthetics)

Prevents flat back

Objective: Prevention of progression of the disease through early application of orthotics.

#### TECHNIQUE

Three C for mobility and high acceptance:

Control – through ongoing correction force

Comfort – through open construction

Cosmetic – through special flat design  
 Trunk movement by hinge

\* Patient is considered from the posterior

Please ask for the dimensions

